

IMPROVE EFFICIENCIES BY REDUCING OPERATING TIME

CORONET ENDOGLIDE™ ULTRATHIN

Tissue is prepared to physician specifications

Graft thickness: 60 - 140 μ

Graft diameter: 7.0 - 8.5mm

Incision Size:

Temporal 4.5mm (scleral tunnel) or

4.9mm (clear corneal) primary incision

Additional Instrumentation Needed:

Placement forceps*

Anterior chamber maintainer



*23g EndoGlide Placement Titanium Forceps (Ref # 53-951)

