

IMPROVE GRAFT ADHERANCE AND ENHANCE SURGICAL CONTROL

MODIFIED JONES TUBE

The Modified Jones Tube Cannula offers several advantages in DMEK surgery, including improved graft adherence and enhanced surgical control. Surgeons can leverage this device to optimize patient outcomes and achieve successful DMEK results.

The tissue has been prepared in an aseptic manner using the modified scuba technique, “F” or “S” stamped (optional), trephinated with a corneal punch, stained with trypan blue and loaded into the glass cannula.

Typical pre-punches are 7.5mm, 7.75mm and 8mm

Approximate Incision Size 3.2mm

Glass cannula is packaged ready to attach to a syringe via luer lock connection.

