

# IMPROVE EFFICIENCIES IN THE OPERATING ROOM AND MINIMIZE TISSUE LOSS


## MICRO JONES TUBE

The tissue has been prepared in an aseptic manner using the modified scuba technique, “F” or “S” stamped (optional), trephinated with a corneal punch, stained with trypan blue and loaded into the glass cannula.

Typical pre-punches are 7.5mm, 7.75mm and 8mm

Approximate Incision Size 2.0-2.2

Glass cannula is packaged ready to attach to a syringe via luer lock connection.



Narrow insertion tip  
allows for a smaller  
incision size

Wide body to protect  
endothelial cells

