

SINGLE LAYER AMNION MEMBRANE FOR OPHTHALMIC USE

DCI AMNIOTIC MEMBRANE GRAFTS

Minimally processed to preserve the native structure of the tissue, dehydrated, and terminally sterilized with gamma irradiation with sterility assurance level (SAL) 10^{-6} to ensure recipient safety. The amnion layer is gently separated from the chorion layer producing a thin, transparent, elastic membrane that acts as a protective barrier.

Maternal donors are screened per FDA and AATB guidelines for a panel of infectious diseases at an FDA registered and CLIA Certified laboratory and against established medical-social risk factors.

ATTRIBUTES

Available as 9mm, 12mm, and 15mm disc.

5-year shelf-life stored at ambient temperature.

Validated industry standard Tyvek and Clearfoil packaging.



Sierra
Donor Services
Eye Bank

DCI
Donor Services
Tissue Bank