

ETHANOL-PRESERVED TISSUE FOR OCULOPLASTIC, GLAUCOMA, AND ANTERIOR SEGMENT SURGERIES

SCLERA

Scleral grafts are preserved in 95% ethanol and stored at room temperature, with a two year expiration time frame.

Ordering specifications:

Scleral grafts are available in 1/4, 1/2, or whole cup sizes.

Custom sizes can be made available upon request.

Guidelines for use:

The sclera must be washed, rinsed and soaked in sterile Balanced Salt Solution (BSS) or sterile Normal Saline (NS) for a minimum of 10 minutes prior to use.

- Re-hydrate the sclera in BSS or NS (60 to 100 cc) for a minimum of 5 minutes.
- Transfer sclera and soak in BSS or NS (60 to 100 cc) for an additional 5 minutes.
- Antibiotic may be added to the BSS or NS if desired.

